



Government of Assam  
**OFFICE OF THE PRINCIPAL-CUM-CHIEF SUPERINTENDENT:: GAUHATI  
MEDICAL COLLEGE & HOSPITAL:: GUWAHATI-32**

No. GMC/MICRO/JE/2021/ 597

Dated: 29/9/2021

**CORRIGENDUM**

A partial modification of the bidding documents of tender vide No. GMC/MICRO/JE/2021/545 Dtd.25.09.2021 for supply of various equipments under the DHR-ICMR funded project entitled, "Understanding the cross talk between the host and the pathogen: A way to identify the novel biomarker for the diagnosis of Japanese Encephalitis virus infection" has been made as per minutes of the pre-bid meeting held on 27.09.2021 in the Department of Microbiology, Gauhati Medical College, and may be downloaded from the website <http://gmchassam.gov.in>.

**Others will remain same.**

  
Dr. Ajanta Sharma, MD, PhD  
Professor  
Department of Microbiology  
Gauhati Medical College

  
Principal-cum-chief Superintendent  
Gauhati Medical College  
Guwahati - 32  
Principal cum-chief Superintendent  
Gauhati Medical College & Hospital  
Guwahati-32

Memo No: GMC/MICRO/JE/2021/ 597

Dated: 29/9/2021

Copy forwarded to:

1. Gauhati Medical College Notice Board
2. Gauhati Medical College & Hospital Notice Board
3. Notice board of Superspeciality Hospital, GMCH
4. Notice board of State Cancer Institute, GMCH
5. Gauhati Medical College website [www.gmchassam.gov.in](http://www.gmchassam.gov.in)
6. Notice board of Office of the Director Medical Education (DME), Govt. of Assam
7. State Public Procurement Portal
8. Office file

  
Principal-cum-chief Superintendent  
Gauhati Medical College  
Guwahati - 32  
Principal cum-chief Superintendent  
Gauhati Medical College & Hospital  
Guwahati-32



Government of Assam

OFFICE OF THE PRINCIPAL-CUM-CHIEF SUPERINTENDENT:: GAUHATI  
MEDICAL COLLEGE & HOSPITAL:: GUWAHATI-32

No. GMC/MICRO/2021/JE/ 598

Dated: 29/7/2021

CORRIGENDUM FOR WEBSITE PUBLICATION

1. Refrigerated centrifuge

Previous specification	Suggested modification	Suggested by
Point 3. Programmability: 3 direct program buttons, plus 96 additional programs accessible via folder	Program: 3-5 direct program button, plus 50 or more programs accessories via folder	Eppendof India Pvt. Ltd.
Point 4. Digital display showing rpm, RCF, time	Digital display showing RPM/ RCF, time and temperature	Eppendof India Pvt. Ltd.
Point 9. Adaptors for 0.5 ml, 0.2 ml tubes and 15 ml conical tube	Adaptors for 0.5 ml, 0.2 adapter and rotor for 15/50ml	Eppendof India Pvt. Ltd.
Point 12. Biocontainment certified sealed buckets	Biocontainment/ autoclavable certified sealed buckets	Eppendof India Pvt. Ltd.


2. Gel Documentation System


Previous specification	Suggested modification	Suggested by
Point 3. Image resolution > 4 megapixel.	Image resolution >15 megapixel.	Eppendof India Pvt. Ltd.
Point 8. Motorized zoom lens with C-mount, f/1.2, 12-75 mm.	Motorized zoom lens with C-mount, f/0.7 or better, 12-75 mm.	Eppendof India Pvt. Ltd.
Point 11. Minimum sample size accommodated should be 28x36 cm	Minimum sample size (gel viewing area) accommodate should be 26 x 21 cm or better	Eppendof India Pvt. Ltd.

  
Dr. Ajanta Sharma, MD, PhD  
Professor  
Department of Microbiology  
Gauhati Medical College

### 3. Multiplex Bead based Cytokine Assay Platform

Previous specification	Suggested modification	Suggested by
Validation kit should be provided with system for IQ/OQ and monthly validation of Optics, Fluidics, Classifier and Reporter channel for ensuring the instrument is running as per specification.	Validation /Combination kit (Preferably validation kit) should be provided with system for IQ/OR for ensuring that the instrument is running as per specifications. Relevant technical notes/document must be attached stating the use and application of validation kit.	ThermoFisher Scientific
Blank polystyrene and magnetic beads alongwith amine coupling kit must be available with the vendor and should be quoted in optional for <u>developing in-house assays</u> . Please attach complete price list in USD and Indian Rupees for these beads and coupling kits.	Blank polystyrene and magnetic beads or other beads coupling kits must be available with the vendor and should be quoted in optional for <u>developing in-house assays</u> . Please attach complete price list in USD and Indian Rupees for these beads and coupling kits.	ThermoFisher Scientific
Should have x-plex platform for designing of pre-mixed custom panels.	Should have x-plex/ procarta plex for designing of pre-mixed custom panels	ThermoFisher Scientific
<b>Software</b> - Performs automated instrument Calibration, Validation and maintenance	Performs automated instrument Calibration, Validation or verification and maintenance	ThermoFisher Scientific
-Should have weighted five parameter SPL logistic equation developed by Brendan Scientific, the leader for ELISA analysis.	Should have weighted five parameter SPL logistic equation	ThermoFisher Scientific

  
**Dr. Ajanta Sharma, MD, PhD**  
 Professor  
 Department of Microbiology  
 Gauhati Medical College

  
 (Dr. A.C. Baishya)  
 Principal-cum-chief Superintendent  
 Gauhati Medical College  
 Guwahati - 32  
 Principal cum-chief Superintendent  
 Gauhati Medical College & Hospital  
 Guwahati-32